



EUROLINE myositis profiles

The most comprehensive line blots for differentiated serological myositis diagnostics



Testing for up to 20 myositis-relevant autoantibodies using one blot strip



Ha, Ks and Zo expand the antigen spectrum for anti-synthetase syndrome (ASS) diagnostics



Detection of anti-cN-1A autoantibodies – the only marker for inclusion body myositis (IBM) diagnostics



Automated processing and evaluation using EUROBlotOne and EUROLineScan



Learn more here:

<http://www.euroimmun.com/myositis>

Example: EUROLINE Autoimmune Inflammatory Myopathies 20 Ag (IgG)

