




Serological assays for parasitic infections



- Broad panel of CE-marked test systems
- Sensitive screening assays
- Careful selection of highly suitable antigens
- Standardised incubation conditions for each test method
- Automated processing possible
- Useful addition to direct detection methods (if indicated)

Pathogen	IIFT	ELISA	BLOT
<i>Echinococcus</i> *		■	■
<i>Leishmania donovani</i>	■		
<i>Plasmodium</i>		■	
<i>Schistosoma</i>	■	■	
<i>Strongyloides</i>		■	
<i>Toxocara canis</i>		■	■
<i>Toxoplasma gondii</i>	■	■	
<i>Trypanosoma cruzi</i>		■	

* ELISA: *Echinococcus* spp., BLOT: *Echinococcus* spp.; differentiation between *E. granulosus* and *E. multilocularis*